

EVACUATION SLYDE®

DQE®

Ready now for what's next.™

FAQs

- Q.** How much do the DQE Evacuation Slyde products weigh?
- A.** The Evacuation Slyde weighs just 9 pounds. The Bariatric Evacuation Slyde weighs 11 pounds.
-
- Q.** What material is used in the Evacuation Slyde and Bariatric Evacuation Slyde?
- A.** They are constructed of durable, yet light, corrugated polypropylene, up to 12 millimeters thick.
-
- Q.** How do I store DQE Evacuation Slyde products so I can access them quickly?
- A.** DQE Evacuation Slydes fold flat and can be stacked. Up to 5 Standard Evacuation Slydes can be stored in the wall-mounted DQE Evacuation Slyde Sleeve. Up to 5 of the wider, Bariatric Evacuation Slydes can be stored in the DQE Bariatric Evacuation Slyde Sleeve. Evacuation Slyde products must be stored at 50° – 90° F.
-
- Q.** How many people are required to control the stair descent of an Evacuation Slyde?
- A.** Based on the weight of the evacuee, a minimum of 2 to 7 people can control the stair descent of an evacuee. For more information, please refer to the operations manual.
-
- Q.** Is the DQE Evacuation Slyde reusable?
- A.** The intended use for the Evacuation Slyde is for a single evacuation incident. However, the Evacuation Slyde can be reused during facility training exercises if prior to each use the Evacuation Slyde is checked for damage and functionality. Should the device ever come in contact with fluids, it should not be reused.
-
- Q.** Does DQE offer training for the Evacuation Slyde?
- A.** Instructional videos are available at shop.dqeready.com/slyde-evacuation-sled/#video. To learn more, contact DQE at info@dqeready.com or 800-355-4628.
-
- Q.** What is the rating of the DQE Stairwell Belay System?
- A.** The equipment is rated at 800 pounds. The system uses 25 feet (7.62 meters) of rope, two carabiners, and a belay friction device to assist with controlled descent.
-
- Q.** What is the purpose of the DQE Stairwell Belay Box?
- A.** Installed on the wall at the top of a flight of stairs, the Belay Box provides convenient storage for the Stairwell Belay System. When properly installed, the included wall anchor provides a secure anchor point for the Stairwell Belay System when other anchor options such as hand railings are unavailable.



CONTACT US

DQE Inc.
www.EvacuationSlyde.com
info@dqeready.com
800-355-4628